



## Screen Detection Levels

Screen ID	Drug/Metabolite	Detection Level
COC*	Benzoylcegonine	300 ng/ml
OPI*	Morphine	300 ng/ml
MET*	d-Methamphetamine	500 ng/ml
THC*	11-nor-79-Tetrahydrocannabinol 9-carboxylic acid	50 ng/ml
AMP*	d-Amphetamine	1000 ng/ml
PCP*	Phencyclidine	25 ng/ml
BZO*	Oxazepam	300 ng/ml
BAR*	Secobarbital	300 ng/ml
MTD*	Methadone	300 ng/ml
MDMA*	3,4 - Methylenedioxymethamphetamine	500 ng/ml
OXY*	Oxycodone	100 ng/ml
TCA**	Nortriptyline	1000 ng/ml
BUP***	Buprenorphine	10 ng/ml

\* 11-panel QuickTox®

\*\* Fastect® II

\*\*\* Fastect® II (available Sept 2008)

Reference: QuickTox® and Fastect® II Reference Guide, Branan Medical Corp.

Screen ID	Common Names	Detection Period
COC	Cocaine, Crack Cocaine	24-60 hours
OPI	Codeine, Morphine, Heroin	3 days
MET	Desoxyn, Meth, Crystal Meth, Ice	5 days
THC	Marijuana	1-15 days
AMP	Biphphetamine, Dexedrine, Speed	5 days
PCP	Phencyclidine, PCP, Angel Dust	3 days
BZO	Ativan, Halcion, Librium, Valium, Xanax	7 days
BAR	Amytal, Nembutal, Seconal, Phenobarbital	6 days
MTD	Methadone	14 days
MDMA	Ecstasy, XTC, Clarity	3-5 days
OXY	OxyContin, Percodan, Percocet, Vicodin, Lortab, Lorcet, Darvon, Darvocet	Varies based on frequency of use, metabolic rate, weight, age, and diet.
TCA	Tricyclic Antidepressants	10 days
BUP	Suboxone	11 days